


Polar Map Templates


When regional polar map or scoring values are exported from Corridor4DM (XML, Reports, Stats), the regional values are ordered based on the overlay as described below.

VASCULAR OVERLAY




Sectors: 1-3 as shown, 4=Total

VASC(NG1) OVERLAY




Sectors: 1-3 as shown, 4=Total

VASC(NG2) OVERLAY




Sectors: 1-3 as shown, 4=Total

5-ZONE OVERLAY




Sectors: 1-5 as shown, 6=Total

9-ZONE OVERLAY




Sectors: 1-9 as shown, 10=Total

13-ZONE OVERLAY




Sectors: 1-13 as shown, 14=Total

17-ZONE OVERLAY




Sectors: 1-17 as shown, 18=Total

19-ZONE OVERLAY



Sectors: 1-19 as shown, 20=Total

20-ZONE OVERLAY



Sectors: 1-20 as shown, 21=Total